


# CANezy

## SPERM SEPARATION DEVICE

CANezy is used to separate motile sperm from semen sample in the easiest way for assisted reproductive technology (ART) procedures like Intrauterine insemination (IUI), In-vitro fertilization (IVF) and Intracytoplasmic sperm injection (ICSI).

### How it works?



Media/ Outlet Port

Sample Port



Volume Adjustment Dial For Sample 0.5ml to 3.0ml



Take desired volume of liquified Semen Sample and inject slowly from Sample port



Take 0.7ml of Sperm Washing/ Flushing Media and Inject through Media/Outlet port slowly

Keep device horizontally & incubate it at 37°C for 15 to 30 minutes.



Aspirate and use for further ART procedure



### KEY FEATURES

- Single device- variable volume (0.5ml to 3.0ml)
- Separation of motile, morphologically normal sperm with low DFI
- Close device- reduces contamination & spills
- Works with fresh & cryopreserved sample
- Volume makeup is only required when semen volume is less than 0.5ml

## Single Device- Variable Volume



**CANDORE BIOSCIENCE**


5, Prime Industrial Park-2, Santej-Vadsar Road,  
Santej-382721, Gujarat, India

[www.candore.in](http://www.candore.in) | [info@candore.in](mailto:info@candore.in) |

+91 - 98240 57910



# FREQUENTLY ASKED QUESTIONS



## **What are the advantages of using CANezy compared to similar devices?**

Variable volume in single device (0.5-3.0 ml), Close device- Reduces evaporation of the sperm preparation medium  
No risk of spillage and contamination

## **On what principle CANezy device works?**

CANezy- Sperm separation device is a close system integrated with microporous filter membrane that allows the sperm to swim through the filter and the separated sperms are then used for ART procedures.

## **Is it compulsory to incubate the sample for 30 minutes?**

No. The required incubation time depends on the concentration of progressively motile sperm in the semen sample. Samples with higher motility generally require a shorter incubation period.

## **What is the sperm concentration limit for using CANezy?**

Minimum sperm concentration of 2 million/mL is required; there is no specified upper limit.

## **What should be the ratio of sperm sample to the media?**

CANezy is the device that allows adjustment of the semen sample volume from 0.5 mL to 3 mL, while maintaining a fixed 0.7 mL volume of separation medium.

## **What sample type can be used for separation of sperm?**

Fresh (raw) or frozen-thawed sample.

## **Is it necessary to wash frozen thaw sample before separation with CANezy?**

Dilute the sample with media, mix gently and use for the separation. It is recommended to check the sample for progressive motile sperms.

## **Does CANezy works efficiently for viscous sample?**

Elevated viscosity may compromise efficient sperm separation. Dilute the sample with media and try to reduce the viscosity without excessive manipulation.

## **What if the incubation time exceeds more than recommended time?**

If the incubation period exceeds, then there are changes of getting sluggish motile sperms and non-motile sperm. When there is a delay in using the separated sperm, remove the separated sperm from the device rather than incubating for longer time and store in appropriate sterile tube.

## **Incubation condition requirement for CANezy device?**

If incubating in HEPES and/or MOPS buffered media then incubate in non-gassed incubator and if using Bicarbonate buffered media use CO<sub>2</sub> incubator.

## **Are there any chances of getting abnormal morphology sperm after incubation?**

There is a possibility that sperms having round head and small acrosome, yet retaining good motility, may pass through the filter membrane during the incubation period.

## **Are there any training or special skills required for using CANezy?**

CANezy is simple to operate and requires no prior training or specialized skills. The only prerequisite is maintaining aseptic conditions throughout the procedure.