





CANros

REACTIVE OXYGEN SPECIES DETECTION KIT



- Measure free radicles present in raw semen sample
- NBT dye based assay
- Results are observed by nacked eyes
- Color Intensity of the vial depends on the concentration of the free radicles present in the semen sample



Interpretation	Colour Code
Normal (White/Light Pink)	
Low level of ROS (Light Purple)	
Moderate level of ROS (Purple)	
High level of ROS (Dark Purple)	



CANDORE BIOSCIENCE

5, Prime Industrial Park-2, Santej-Vadsar Road, Santej- 383721, Gujarat, INDIA

+91 98240 57910

info@candore.in

www.candore.com