



# PRP SEPARATION KIT

“Engineered for clinicians who expect more.”

- Close system- Minimized risk of contamination
- 4 to 5 fold platelet concentration in the isolate
- One step process- only 6 minutes
- Kit contains Anticoagulant & Activation solution
- Easy to use
- No specialized centrifuge is required

CANDORE BIOSCIENCE

5, Prime Industrial Park-2, Santej- Vadsar Road, Santej- 383721, Gujarat INDIA

+91 98240 57910

info@candore.in

www.candore.in

